


Portable Urine Analyzer

CYBOW™ R-50S





Product Features

- Special urine test such as Microalbuminuria test, Hypercalciuria test etc. are available besides general urine tests using with DFI's urinalysis strips
- Simple operation by Touch LCD
- Light and compact (portable) size
- Automatic display of test results and relevant information on LCD screen
- Result can be obtained after 60 sec.

Specifications

- Measurement : Reflectance Photometer
- Size(Dimension) : 188 (W) x 74 (D) x 77 (H) mm
- Net Weight: 0.40 Kg
- Power : DC 12V / 0.33A or battery 1.5V x 8EA
- Throughput : 45 (general mode) ~ 120 (speed mode) tests per hour
- Interface : USB
- Weight : 430g
- Memory: 2000 results
- Display: 4.3 inch TFT LCD
- Software : user program PC s/w is available
- Result printing : via PC s/w or External printer (optional)



ComboStik[®]

Urine Reagent Strips
Urine Chemistry Analyzer

